

GOLPH2 Rabbit mAb

Catalog No: #58875

Package Size: #58875-1 50ul #58875-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GOLPH2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	GOLPH2 Antibody detects endogenous levels of GOLPH2
Immunogen Description	A synthesized peptide derived from human GOLPH2
Other Names	GOLPH2; GP73; GOLM1; Golgi protein 73kD; C9orf155;
Accession No.	Uniprot:Q8NBJ4
Uniprot	Q8NBJ4
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of GOLPH2 expression in LnCaP cell lysate.

Product Description

Unknown. Cellular response protein to viral infection.

Background

Unknown. Cellular response protein to viral infection.

Note: This product is for in vitro research use only