

## IRF7 Rabbit mAb

Catalog No: #58927

Package Size: #58927-1 50ul #58927-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IRF7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	IRF7 Antibody detects endogenous levels of IRF7
Immunogen Description	A synthesized peptide derived from human IRF7
Other Names	IRF7; Interferon regulatory factor 7; IRF-7; IRF7A; IRF-7H;
Accession No.	Uniprot:Q92985
Uniprot	Q92985
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of IRF7 expression in Jurkat cell lysate.



Immunohistochemical analysis of paraffin-embedded human pancreas, using IRF7 Antibody.

## Product Description

Binds to the Q promoter (Qp) of EBV nuclear antigen 1 a (EBNA1) and may play a role in the regulation of EBV latency. Can activate distinct gene expression programs in macrophages and regulate the anti-tumor properties of primary macrophages.

## Background

Binds to the Q promoter (Qp) of EBV nuclear antigen 1 a (EBNA1) and may play a role in the regulation of EBV latency. Can activate distinct gene expression programs in macrophages and regulate the anti-tumor properties of primary macrophages.

Note: This product is for in vitro research use only