

## IRAK3 Rabbit mAb

Catalog No: #58947

Package Size: #58947-1 50ul #58947-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	IRAK3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	IRAK3 Antibody detects endogenous levels of IRAK3
Immunogen Description	A synthesized peptide derived from human IRAK3
Other Names	ASRT5; FLJ13601; IL-1 receptor-associated kinase M; Interleukin-1 receptor-associated kinase 3; interleukin-1 receptor-associated kinase M; IRAK-3; IRAK-M; IRAK3; IRAKM;
Accession No.	Uniprot:Q9Y616
Uniprot	Q9Y616
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details


WB 1:500~1:2000 IP 1:10

## Images



Western blot analysis of SGK1 expression in A431 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

IRAK3 a TKL kinase of the IRAK family. The interleukin-1 receptor-associated kinases are important mediators in the signal transduction of Toll-like receptor and IL1R family members, collectively referred to as TIRs.

## Background

IRAK3 a TKL kinase of the IRAK family. The interleukin-1 receptor-associated kinases are important mediators in the signal transduction of Toll-like receptor and IL1R family members, collectively referred to as TIRs.

Note: This product is for in vitro research use only