

## CXCR7 Rabbit mAb

Catalog No: #58958

Package Size: #58958-1 50ul #58958-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CXCR7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	CXCR7 Antibody detects endogenous levels of CXCR7
Immunogen Description	A synthesized peptide derived from human CXCR7
Other Names	CMKOR1; RDC1; CXCR7; GPR159; GPRN1; ACKR3;
Accession No.	Uniprot:P25106
Uniprot	P25106
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 FC 1:50

## Images



Western blot analysis of CXCR7 expression in K562 cell lysate.

## Product Description

Plays a regulatory role in CXCR4-mediated activation of cell surface integrins by CXCL12. Required for heart valve development. Acts as coreceptor with CXCR4 for a restricted number of HIV isolates.

## Background

---

Plays a regulatory role in CXCR4-mediated activation of cell surface integrins by CXCL12. Required for heart valve development. Acts as coreceptor with CXCR4 for a restricted number of HIV isolates.

---

Note: This product is for in vitro research use only