

## AGTR2 Rabbit mAb

Catalog No: #58986

Package Size: #58986-1 50ul #58986-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	AGTR2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	AGTR2 Antibody detects endogenous levels of AGTR2
Immunogen Description	A synthesized peptide derived from human AGTR2
Other Names	AGTR2; Angiotensin ii receptor type 2; Angiotensin II type-2 receptor; AT2 receptor; AT2; At2r; MRX88; Angiotensin ii at2; Angiotensin receptor 2; Type-2 angiotensin II receptor; ATGR2;
Accession No.	Uniprot:P50052
Uniprot	P50052
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:2000 IP 1:50

## Images



Western blot analysis of AGTR2 expression in HepG2 cell lysate.

## Product Description

Receptor for angiotensin II. Cooperates with MTUS1 to inhibit ERK2 activation and cell proliferation.

## Background

---

Receptor for angiotensin II. Cooperates with MTUS1 to inhibit ERK2 activation and cell proliferation.

---

Note: This product is for in vitro research use only