

Caspase 5 Rabbit mAb

Catalog No: #59005

Package Size: #59005-1 50ul #59005-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Caspase 5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	Caspase 5 Antibody detects endogenous levels of Caspase 5
Immunogen Description	A synthesized peptide derived from human Caspase 5
Other Names	CASP5; Apoptosis related cysteine protease; CASP-5; Caspase-5 subunit p10; Protease ICH-3; TY protease;
Accession No.	Uniprot:P51878
Uniprot	P51878
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:1000~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of Caspase 5 expression in HeLa cell lysate.



Immunohistochemical analysis of paraffin-embedded human breast cancer, using Caspase 5 Antibody.

Product Description

Caspase-5 (Ich-3/ICErelIII/TY) is a member of the caspase family of cysteine proteases that play a key role in the execution of apoptosis and activation of inflammatory cytokines. It is a mediator of programmed cell death (apoptosis).

Background

Caspase-5 (Ich-3/ICErelIII/TY) is a member of the caspase family of cysteine proteases that play a key role in the execution of apoptosis and activation of inflammatory cytokines. It is a mediator of programmed cell death (apoptosis).

Note: This product is for in vitro research use only