

Histone H3 (di methyl K9) Rabbit mAb

Catalog No: #59023

Package Size: #59023-1 50ul #59023-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Histone H3 (di methyl K9) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF ChIP
Species Reactivity	Human Mouse Rat
Specificity	Histone H3 (di methyl K9) Antibody detects endogenous levels of Histone H3 (di methyl K9)
Immunogen Description	A synthesized peptide derived from human Histone H3
Other Names	Histone H3.1; Histone H3; HIST1H3A; H3K9me2;
Accession No.	Uniprot:P68431
Uniprot	P68431
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 ICC/IF 1:50~1:200 ChIP

Images



Western blot analysis of Histone H3 (di methyl K9) expression in HeLa cell lysate.

Product Description

Belongs to the histone H3 family. Play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

Background

Belongs to the histone H3 family. Play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.

Note: This product is for in vitro research use only