

INPP4B Rabbit mAb

Catalog No: #59032

Package Size: #59032-1 50ul #59032-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	INPP4B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	INPP4B Antibody detects endogenous levels of INPP4B
Immunogen Description	A synthesized peptide derived from human INPP4B
Other Names	4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa; Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3,4 bisphosphate 4 phosphatase; Type II inositol-3;
Accession No.	Uniprot:O15327
Uniprot	O15327
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:2000 IHC 1:50~1:200 FC 1:50

Images



Western blot analysis of INPP4B expression in MCF-7 cell lysate.

Product Description

Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

Background

Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

Note: This product is for in vitro research use only