

TWEAKR Rabbit mAb

Catalog No: #59056

Package Size: #59056-1 50ul #59056-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TWEAKR Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	TWEAKR Antibody detects endogenous levels of TWEAKR
Immunogen Description	A synthesized peptide derived from human TWEAKR
Other Names	CD266 antigen; FGF inducible 14; Fibroblast growth factor inducible immediate early response protein 14; FN14; TNFRSF 12A; Tweak receptor; TweakR;
Accession No.	Uniprot:Q9NP84
Uniprot	Q9NP84
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of TWEAKR expression in HUVEC cell lysate.

Product Description

Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins.

Background

Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins.

Note: This product is for in vitro research use only