

Cytokeratin 9 Rabbit mAb

Catalog No: #59067

Package Size: #59067-1 50ul #59067-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Cytokeratin 9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	Cytokeratin 9 Antibody detects endogenous levels of Cytokeratin 9
Immunogen Description	A synthesized peptide derived from human Cytokeratin 9
Other Names	CK 9; Cytokeratin 9; EPPK; K9; Keratin 9; Keratin; Keratin type I cytoskeletal 9; Keratin9; KRT 9; type I cytoskeletal 9;
Accession No.	Uniprot:P35527
Uniprot	P35527
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of Cytokeratin 9 expression in HeLa cell lysate.

Product Description

May serve an important special function either in the mature palmar and plantar skin tissue or in the morphogenetic program of the formation of these tissues. Plays a role in keratin filament assembly.

Background

May serve an important special function either in the mature palmar and plantar skin tissue or in the morphogenetic program of the formation of these tissues. Plays a role in keratin filament assembly.

Note: This product is for in vitro research use only