

## SOX10 Rabbit mAb

Catalog No: #59083

Package Size: #59083-1 50ul #59083-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	SOX10 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	SOX10 Antibody detects endogenous levels of SOX10
Immunogen Description	A synthesized peptide derived from human SOX10
Other Names	SOX10; DOM; Transcription factor SOX-10; WS4; PCWH; WS2E; SRY-related HMG-box gene 10; WS4C;
Accession No.	Uniprot:P56693
Uniprot	P56693
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

IHC 1:50~1:200

## Images



Immunohistochemical analysis of paraffin-embedded human breast, using SOX10 Antibody.

## Product Description

SOX10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and the determination of cell fate. Sox10 is an important regulator of neural crest and peripheral nervous system development.

## Background

---

SOX10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and the determination of cell fate. Sox10 is an important regulator of neural crest and peripheral nervous system development.

---

Note: This product is for in vitro research use only