

## TRK fused gene Rabbit mAb

Catalog No: #59132

Package Size: #59132-1 50ul #59132-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TRK fused gene Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	TRK fused gene Antibody detects endogenous levels of TRK fused gene
Immunogen Description	A synthesized peptide derived from human TRK fused gene
Other Names	HMSNP; Protein TFG; SPG57; TF6; TFG; TRK fused; TRK fused gene protein; TRKT3; TRKT3 oncogene;
Accession No.	Uniprot:Q92734
Uniprot	Q92734
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:5000~1:20000 IHC 1:50~1:200 ICC/IF 1:100~1:500 IP 1:50

## Images



Western blot analysis of TRK fused gene expression in HeLa cell lysate.

## Product Description

Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

## Background

Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

---

Note: This product is for in vitro research use only