

## PAX9 Rabbit mAb

Catalog No: #59134

Package Size: #59134-1 50ul #59134-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PAX9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	PAX9 Antibody detects endogenous levels of PAX9
Immunogen Description	A synthesized peptide derived from human PAX9
Other Names	Paired box protein Pax-9; PAX9;
Accession No.	Uniprot:P55771
Uniprot	P55771
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:2000

## Images



Western blot analysis of PAX9 expression in Jurkat cell lysate.

## Product Description

Paired box (PAX) proteins are a family of transcription factors that play important and diverse roles in animal development. Transcription factor required for normal development of thymus, parathyroid glands, ultimobranchial bodies, teeth, skeletal elements of skull and larynx as well as distal limbs.

## Background

---

Paired box (PAX) proteins are a family of transcription factors that play important and diverse roles in animal development. Transcription factor required for normal development of thymus, parathyroid glands, ultimobranchial bodies, teeth, skeletal elements of skull and larynx as well as distal limbs.

---

Note: This product is for in vitro research use only