

PAX8 Rabbit mAb

Catalog No: #59142

Package Size: #59142-1 50ul #59142-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | PAX8 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | IHC ICC/IF |
| Species Reactivity | Human |
| Specificity | PAX8 Antibody detects endogenous levels of PAX8 |
| Immunogen Description | A synthesized peptide derived from human PAX8 |
| Other Names | Paired box 8; Paired box gene 8; paired box homeotic gene 8; Paired box protein Pax 8; Paired domain gene 8; PAX 8; |
| Accession No. | Uniprot:Q06710 |
| Uniprot | Q06710 |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

IHC 1:50~1:200 ICC/IF 1:50~1:200

Product Description

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

Background

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

Note: This product is for in vitro research use only