

MGMT Rabbit mAb

Catalog No: #59191

Package Size: #59191-1 50ul #59191-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MGMT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	MGMT Antibody detects endogenous levels of total MGMT
Immunogen Description	A synthesized peptide derived from human MGMT
Other Names	MGMT;
Accession No.	Uniprot:P16455
Uniprot	P16455
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of MGMT expression in HeLa cell lysate.

Product Description

MGMT (O-6-methylguanine-DNA methyltransferase) is a DNA repair enzyme that participates in a suicide reaction that specifically removes methyl or alkyl groups from the O(6) position of guanine, restoring guanine to its normal form without causing DNA breaks. MGMT protects cells from alkylating toxins, and is an important factor in drug resistance to alkylating therapeutic agents.

Background

MGMT (O-6-methylguanine-DNA methyltransferase) is a DNA repair enzyme that participates in a suicide reaction that specifically removes methyl or alkyl groups from the O(6) position of guanine, restoring guanine to its normal form without causing DNA breaks. MGMT protects cells from alkylating toxins, and is an important factor in drug resistance to alkylating therapeutic agents.

Note: This product is for in vitro research use only