

HGF Antibody HRP Conjugated

Catalog No: #C00675H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HGF Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 210-260 723 derived from human HGF alpha
Conjugates	HRP
Target Name	HGF
Other Names	SF; HGFB; HPTA; F-TCF; DFNB39; Hepatocyte growth factor; Hepatopoietin-A; Scatter factor; HGF
Accession No.	Swiss-Prot#P14210NCBI Gene ID3082
Uniprot	P14210
GeneID	3082;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Potent mitogen for mature parenchymal hepatocyte cells, seems to be a hepatotrophic factor, and acts as a growth factor for a broad spectrum of tissues and cell types. Activating ligand for the receptor tyrosine kinase MET by binding to it and promoting its dimerization.

Note: This product is for in vitro research use only