

Collagen X Rabbit mAb

Catalog No: #59247

Package Size: #59247-1 50ul #59247-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Collagen X Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Collagen X Antibody detects endogenous levels of total Collagen X
Immunogen Description	A synthesized peptide derived from human Collagen X
Other Names	CollagenX; COL10A1; COAA1; Collagen type X alpha 1;
Accession No.	Uniprot:Q03692
Uniprot	Q03692
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Collagen X expression in human fetal skin lysate.

Product Description

Type X collagen is a product of hypertrophic chondrocytes and has been localized to presumptive mineralization zones of hyaline cartilage.

Background

Type X collagen is a product of hypertrophic chondrocytes and has been localized to presumptive mineralization zones of hyaline cartilage.

Note: This product is for in vitro research use only