

ACTN2 Rabbit mAb

Catalog No: #59256

Package Size: #59256-1 50ul #59256-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACTN2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	ACTN2 Antibody detects endogenous levels of total ACTN2
Immunogen Description	A synthesized peptide derived from human ACTN2
Other Names	Actin binding protein; Actinin alpha 2; ACTN 2; ACTN2; Alpha actinin 2;
Accession No.	Uniprot:P35609
Uniprot	P35609
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of ACTN2 expression in HepG2 cell lysate.

Product Description

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Background

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Note: This product is for in vitro research use only