

NEFH Rabbit mAb

Catalog No: #59274

Package Size: #59274-1 50ul #59274-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NEFH Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	NEFH Antibody detects endogenous levels of total NEFH
Immunogen Description	A synthesized peptide derived from human NEFH
Other Names	NF-H; 200 kDa neurofilament protein; Neurofilament triplet H protein; KIAA0845; NFH; NEFH;
Accession No.	Uniprot:P12036
Uniprot	P12036
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500

Images



Western blot analysis of NEFH expression in human brain lysate.

Product Description

Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber. NF-H has an important function in mature axons that is not subserved by the two smaller NF proteins.

Background

Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber. NF-H has an important function in mature axons that is not subserved by the two smaller NF proteins.

Note: This product is for in vitro research use only