

## TBR1 Rabbit mAb

Catalog No: #59277

Package Size: #59277-1 50ul #59277-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TBR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	TBR1 Antibody detects endogenous levels of total TBR1
Immunogen Description	A synthesized peptide derived from human TBR1
Other Names	T box brain 1; T-brain-1; TBR 1; Tbr1; TES 56;
Accession No.	Uniprot:Q16650
Uniprot	Q16650
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of TBR1 expression in human fetal brain lysate.

## Product Description

Probable transcriptional regulator involved in developmental processes. Required for normal brain development.

## Background

Probable transcriptional regulator involved in developmental processes. Required for normal brain development.

---

Note: This product is for in vitro research use only