

## Parvalbumin Rabbit mAb

Catalog No: #59286

Package Size: #59286-1 50ul #59286-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Parvalbumin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Parvalbumin Antibody detects endogenous levels of total Parvalbumin
Immunogen Description	A synthesized peptide derived from human Parvalbumin
Other Names	Parvalbumin alpha; PVALB;
Accession No.	Uniprot:P20472
Uniprot	P20472
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

## Images



Western blot analysis of Parvalbumin expression in human cerebellum lysate.

## Product Description

In muscle, parvalbumin is thought to be involved in relaxation after contraction. It binds two calcium ions.

## Background

In muscle, parvalbumin is thought to be involved in relaxation after contraction. It binds two calcium ions.

---

Note: This product is for in vitro research use only