

## Prealbumin Rabbit mAb

Catalog No: #59292

Package Size: #59292-1 50ul #59292-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Prealbumin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Prealbumin Antibody detects endogenous levels of total Prealbumin
Immunogen Description	A synthesized peptide derived from human Prealbumin
Other Names	Amyloid polyneuropathy; Amyloidosis I; ATTR; CTS; PALB; Prealbumin amyloidosis type I; Prealbumin;
Accession No.	Uniprot:P02766
Uniprot	P02766
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

## Images



Western blot analysis of Prealbumin gamma expression in human fetal kidney lysate.

## Product Description

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

## Background

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

---

Note: This product is for in vitro research use only