

PROX1 Rabbit mAb

Catalog No: #59294

Package Size: #59294-1 50ul #59294-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PROX1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	PROX1 Antibody detects endogenous levels of total PROX1
Immunogen Description	A synthesized peptide derived from human PROX1
Other Names	Homeobox prospero like protein PROX1; Prospero homeobox protein 1; prospero-related homeobox gene 1; PROX1;
Accession No.	Uniprot:Q92786
Uniprot	Q92786
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of PROX1 expression in HepG2 cell lysate.

Product Description

Plays a critical role in embryonic development and functions as a key regulatory protein in neurogenesis and the development of the heart, eye lens, liver, pancreas and the lymphatic system. Involved in the regulation of the circadian rhythm.

Background

Plays a critical role in embryonic development and functions as a key regulatory protein in neurogenesis and the development of the heart, eye lens, liver, pancreas and the lymphatic system. Involved in the regulation of the circadian rhythm.

Note: This product is for in vitro research use only