

## AGTR1 Rabbit mAb

Catalog No: #59297

Package Size: #59297-1 50ul #59297-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	AGTR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	AGTR1 Antibody detects endogenous levels of total AGTR1
Immunogen Description	A synthesized peptide derived from human AGTR1
Other Names	AGTR1; AG2S; AGTR1A; AGTR1B; AT1; AT1AR; AT1B; AT1BR; AT2R1; AT2R1A; AT2R1B; HAT1R;
Accession No.	Uniprot:P30556
Uniprot	P30556
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of AGTR1 expression in (1) HeLa cell lysate; (2) PC-12 cell lysate.

## Product Description

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

## Background

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

---

Note: This product is for in vitro research use only