

## TGN46 Rabbit mAb

Catalog No: #59302

Package Size: #59302-1 50ul #59302-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TGN46 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	TGN46 Antibody detects endogenous levels of total TGN46
Immunogen Description	A synthesized peptide derived from human TGN46
Other Names	TGN38; TGN46; TGN48; TGN51; TGOLN2; TTGN2;
Accession No.	Uniprot:O43493
Uniprot	O43493
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

## Images



Western blot analysis of TGN46 expression in HeLa cell lysate.

## Product Description

May be involved in regulating membrane traffic to and from trans-Golgi network.

## Background

May be involved in regulating membrane traffic to and from trans-Golgi network.

---

Note: This product is for in vitro research use only