

GC / VDBP Rabbit mAb

Catalog No: #59313

Package Size: #59313-1 50ul #59313-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GC / VDBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	GC / VDBP Antibody detects endogenous levels of total GC / VDBP
Immunogen Description	A synthesized peptide derived from human GC / VDBP
Other Names	DBP; DBP/GC; GRD3; VDBG; VDBP;
Accession No.	Uniprot:P02774
Uniprot	P02774
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:30

Images



Western blot analysis of GC/VDBP expression in Human plasma lysate.

Product Description

DBP associates with membrane-bound immunoglobulin on the surface of B-lymphocytes and with IgG Fc receptor on the membranes of T-lymphocytes.

Background

DBP associates with membrane-bound immunoglobulin on the surface of B-lymphocytes and with IgG Fc receptor on the membranes of T-lymphocytes.

Note: This product is for in vitro research use only