

## Chromogranin A Rabbit mAb

Catalog No: #59314

Package Size: #59314-1 50ul #59314-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Chromogranin A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	Chromogranin A Antibody detects endogenous levels of total Chromogranin A
Immunogen Description	A synthesized peptide derived from human Chromogranin A
Other Names	beta Granin; Chromogranin A; ChrA; SP-I; Vasostatin I; Vasostatin-2; Vasostatin II; Pancreastatin; CHGA;
Accession No.	Uniprot:P10645
Uniprot	P10645
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:20

## Images



Western blot analysis of Chromogranin A expression in SH-SY5Y cell lysate.

## Product Description

Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

## Background

Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

---

Note: This product is for in vitro research use only