

## Urokinase Rabbit mAb

Catalog No: #59337

Package Size: #59337-1 50ul #59337-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Urokinase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Urokinase Antibody detects endogenous levels of total Urokinase
Immunogen Description	A synthesized peptide derived from human Urokinase
Other Names	Urokinase-type plasminogen activator; U-plasminogen activator; Upa; PLAU; ATF uPA;
Accession No.	Uniprot:P00749
Uniprot	P00749
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:1000 IHC 1:50~1:200

## Images



Western blot analysis of Urokinase expression in LnCaP cell lysate.

## Product Description

Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.

## Background

Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.

---

Note: This product is for in vitro research use only