

FGFR3 Rabbit mAb

Catalog No: #59348

Package Size: #59348-1 50ul #59348-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FGFR3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	FGFR3 Antibody detects endogenous levels of total FGFR3
Immunogen Description	A synthesized peptide derived from human FGFR3
Other Names	Fibroblast growth factor receptor 3; CD333; FGFR3; JTK4; MFR 3; SAM 3;
Accession No.	Uniprot:P22607
Uniprot	P22607
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of FGFR3 expression in HeLa cell lysate.

Product Description

Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

Background

Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

Note: This product is for in vitro research use only