

IP10 Rabbit mAb

Catalog No: #59349

Package Size: #59349-1 50ul #59349-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IP10 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IP10 Antibody detects endogenous levels of total IP10
Immunogen Description	A synthesized peptide derived from human IP10
Other Names	C-X-C motif chemokine 10; 10 kDa interferon gamma-induced protein; Gamma-IP10; IP-10; Small-inducible cytokine B10; CXCL10; INP10; SCYB10; MOB1;
Accession No.	Uniprot:P02778
Uniprot	P02778
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of IP10 expression in IP10 recombinant protein lysate.

Product Description

Chemotactic for monocytes and T-lymphocytes. Binds to CXCR3.

Background

Chemotactic for monocytes and T-lymphocytes. Binds to CXCR3.

Note: This product is for in vitro research use only