

Perilipin A Rabbit mAb

Catalog No: #59362

Package Size: #59362-1 50ul #59362-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Perilipin A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse
Specificity	Perilipin A Antibody detects endogenous levels of total Perilipin A
Immunogen Description	A synthesized peptide derived from human Perilipin A
Other Names	Lipid droplet associated protein; PERI; Perilipin-1; PerilipinA; LIN1;
Accession No.	Uniprot:O60240
Uniprot	O60240
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 ICC/IF 1:50~1:200

Images



Western blot analysis of Perilipin A expression in human fetal liver lysate.

Product Description

Modulator of adipocyte lipid metabolism. Coats lipid storage droplets to protect them from breakdown by hormone-sensitive lipase (HSL). Its absence may result in leanness.

Background

Modulator of adipocyte lipid metabolism. Coats lipid storage droplets to protect them from breakdown by hormone-sensitive lipase (HSL). Its absence may result in leanness.

Note: This product is for in vitro research use only