

FGF21 Rabbit mAb

Catalog No: #59363

Package Size: #59363-1 50ul #59363-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FGF21 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	FGF21 Antibody detects endogenous levels of total FGF21
Immunogen Description	A synthesized peptide derived from human FGF21
Other Names	Fgf21; FGFL; Fibroblast growth factor 21; PRO10196; UNQ3115;
Accession No.	Uniprot:Q9NSA1
Uniprot	Q9NSA1
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of FGF21 expression in (1) HeLa cell lysate; (2) Mouse spleen lysate.

Product Description

Stimulates glucose uptake in differentiated adipocytes via the induction of glucose transporter SLC2A1/GLUT1 expression (but not SLC2A4/GLUT4 expression). Activity requires the presence of KLB.

Background

Stimulates glucose uptake in differentiated adipocytes via the induction of glucose transporter SLC2A1/GLUT1 expression (but not SLC2A4/GLUT4 expression). Activity requires the presence of KLB.

Note: This product is for in vitro research use only