

## FABP4 Rabbit mAb

Catalog No: #59365

Package Size: #59365-1 50ul #59365-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FABP4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	FABP4 Antibody detects endogenous levels of total FABP4
Immunogen Description	A synthesized peptide derived from human FABP4
Other Names	422/aP2; adipocyte; Adipocyte lipid binding protein; ALBP; ALBP/Ap2; aP2; FABP; FABP4;
Accession No.	Uniprot:P15090
Uniprot	P15090
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of FABP4 expression in human fetal heart lysate.

## Product Description

Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus.

## Background

---

Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus.

---

Note: This product is for in vitro research use only