

NFAT1 Rabbit mAb

Catalog No: #59375

Package Size: #59375-1 50ul #59375-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NFAT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	26A12
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	NFAT1 Antibody detects endogenous levels of total NFAT1
Immunogen Description	A synthesized peptide derived from human NFAT1
Target Name	NFAT1
Other Names	cytoplasmic 2; NF ATp; NF-ATc2; NF-ATp; NFAT 1; NFAT1; NFAT1-D; NFATc2; NFATp;
Accession No.	Uniprot:Q13469
Uniprot	Q13469
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



All Lanes:NFAT1 Antibody at 1/ 1000 dilution.Lane1:Ramos cell lysate.Proteins at 20ug per laneSecondary: Goat Anti-Rabbit IgG(HRP) at 1/20000 dilutionPredicted band size : 135kDaObserved band size: 135kDa

Product Description

Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2, IL-3, IL-4, TNF-alpha or GM-CSF.

Background

Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2, IL-3, IL-4, TNF-alpha or GM-CSF.

Note: This product is for in vitro research use only