

GAD67 Rabbit mAb

Catalog No: #59378

Package Size: #59378-1 50ul #59378-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GAD67 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	GAD67 Antibody detects endogenous levels of total GAD67
Immunogen Description	A synthesized peptide derived from human GAD67
Other Names	Glutamate decarboxylase 1; 67 kDa glutamic acid decarboxylase; GAD 67; GAD1; GAD; GAD67;
Accession No.	Uniprot:Q99259
Uniprot	Q99259
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:50

Images



Western blot analysis of GAD67 expression in (1) Jurkat cell lysate; (2) C6 cell lysate.

Product Description

Catalyzes the production of GABA.

Background

Catalyzes the production of GABA.

Note: This product is for in vitro research use only