

ADFP Rabbit mAb

Catalog No: #59404

Package Size: #59404-1 50ul #59404-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ADFP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	ADFP Antibody detects endogenous levels of total ADFP
Immunogen Description	A synthesized peptide derived from human ADFP
Other Names	ADFP; Adipophilin; Adipose differentiation related protein; ADRP; Perilipin-2; PLIN2;
Accession No.	Uniprot:Q99541
Uniprot	Q99541
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of ADFP expression in HepG2 cell lysate.

Product Description

May be involved in development and maintenance of adipose tissue.

Background

May be involved in development and maintenance of adipose tissue.

Note: This product is for in vitro research use only