

CD73 Rabbit mAb

Catalog No: #59406

Package Size: #59406-1 50ul #59406-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD73 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CD73 Antibody detects endogenous levels of total CD73
Immunogen Description	A synthesized peptide derived from human CD73
Other Names	5' nucleotidase (CD73); 5' nucleotidase precursor; 5' nucleotidase; ecto; 5NTD; CD73; CD73 antigen; E5NT; Ecto 5' nucleotidase; Ecto-5'-nucleotidase; NT5E; NTE;
Accession No.	Uniprot:P21589
Uniprot	P21589
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of CD73 expression in A375 cell lysate.

Product Description

Hydrolyzes extracellular nucleotides into membrane permeable nucleosides.

Background

Hydrolyzes extracellular nucleotides into membrane permeable nucleosides.

Note: This product is for in vitro research use only