

## TTR Rabbit mAb

Catalog No: #59409

Package Size: #59409-1 50ul #59409-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TTR Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	TTR Antibody detects endogenous levels of total TTR
Immunogen Description	A synthesized peptide derived from human TTR
Other Names	Amyloid polyneuropathy; Amyloidosis I; ATTR; Carpal tunnel syndrome 1; CTS; CTS1; HsT2651; PALB; Prealbumin amyloidosis type I; Prealbumin;
Accession No.	Uniprot:P02766
Uniprot	P02766
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of TTR expression in human serum lysate.

## Product Description

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

## Background

---

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

---

Note: This product is for in vitro research use only