

GDNF Rabbit mAb

Catalog No: #59438

Package Size: #59438-1 50ul #59438-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GDNF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	GDNF Antibody detects endogenous levels of total GDNF
Immunogen Description	A synthesized peptide derived from human GDNF
Other Names	ATF1; ATF2; gdnf; HFB1 GDNF; hGDNF; HSCR3;
Accession No.	Uniprot:P39905
Uniprot	P39905
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of GDNF expression in HepG2 cell lysate.

Product Description

Neurotrophic factor that enhances survival and morphological differentiation of dopaminergic neurons and increases their high-affinity dopamine uptake.

Background

Neurotrophic factor that enhances survival and morphological differentiation of dopaminergic neurons and increases their high-affinity dopamine uptake.

Note: This product is for in vitro research use only