

## Involucrin Rabbit mAb

Catalog No: #59444

Package Size: #59444-1 50ul #59444-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Involucrin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Rat
Specificity	Involucrin Antibody detects endogenous levels of total Involucrin
Immunogen Description	A synthesized peptide derived from human Involucrin
Other Names	Involucrin; IVL;
Accession No.	Uniprot:P07476
Uniprot	P07476
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 FC 1:10

## Images



Western blot analysis of Involucrin expression in A431 cell lysate.

## Product Description

Part of the insoluble cornified cell envelope (CE) of stratified squamous epithelia.

## Background

Part of the insoluble cornified cell envelope (CE) of stratified squamous epithelia.

---

Note: This product is for in vitro research use only