

TGM2 Rabbit mAb

Catalog No: #59449

Package Size: #59449-1 50ul #59449-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TGM2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse
Specificity	TGM2 Antibody detects endogenous levels of total TGM2
Immunogen Description	A synthesized peptide derived from human TGM2
Other Names	TGM2; Tissue transglutaminase; Transglutaminase C; TG(C); TGC; TGase C; ransglutaminase H; TGase H; Transglutaminase-2; TGase-2;
Accession No.	Uniprot:P21980
Uniprot	P21980
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of TGM2 expression in HUVEC cell lysate.

Product Description

Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

Background

Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

Note: This product is for in vitro research use only