

DR5 Rabbit mAb

Catalog No: #59451

Package Size: #59451-1 50ul #59451-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DR5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	DR5 Antibody detects endogenous levels of total DR5
Immunogen Description	A synthesized peptide derived from human DR5
Other Names	Death receptor 5; TRAIL receptor 2; TRAIL-R2; CD262; TNFRSF10B; DR5; KILLER; TRAILR2; TRICK2; ZTNFR9;
Accession No.	Uniprot:O14763
Uniprot	O14763
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IP 1:30 FC 1:100

Images



Western blot analysis of DR5 expression in HeLa cell lysate.

Product Description

Plays an important role in the regulation of apoptosis in various physiological systems. Promotes the activation of NF-kappa-B. Essential for ER stress-induced apoptosis.

Background

Plays an important role in the regulation of apoptosis in various physiological systems. Promotes the activation of NF-kappa-B. Essential for ER stress-induced apoptosis.

Note: This product is for in vitro research use only