

GLUT3 Rabbit mAb

Catalog No: #59462

Package Size: #59462-1 50ul #59462-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GLUT3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	GLUT3 Antibody detects endogenous levels of total GLUT3
Immunogen Description	A synthesized peptide derived from human GLUT3
Other Names	brain; Glucose transporter type 3; GLUT3; Slc2a3;
Accession No.	Uniprot:P11169
Uniprot	P11169
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of GLUT3 expression in (1) HepG2 cell lysate; (2) NIH/3T3 cell lysate.

Product Description

Facilitative glucose transporter. Probably a neuronal glucose transporter.

Background

Facilitative glucose transporter. Probably a neuronal glucose transporter.

Note: This product is for in vitro research use only