

Apg12(Atg12) Rabbit mAb

Catalog No: #59469

Package Size: #59469-1 50ul #59469-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Apg12(Atg12) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	Apg12(Atg12) Antibody detects endogenous levels of total Apg12(Atg12)
Immunogen Description	A synthesized peptide derived from human Apg12(Atg12)
Other Names	APG12-like; APG12L; ATG12; Autophagy 12; Autophagy-related protein 12; FBR93; HAPG12; Ubiquitin-like protein ATG12;
Accession No.	Uniprot:O94817
Uniprot	O94817
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:100

Images



Western blot analysis of Apg12(Atg12) expression in 293 cell lysate.

Product Description

Ubiquitin-like protein required for autophagy. Conjugated to ATG3 and ATG5.

Background

Ubiquitin-like protein required for autophagy. Conjugated to ATG3 and ATG5.

Note: This product is for in vitro research use only