

## Periostin Rabbit mAb

Catalog No: #59472

Package Size: #59472-1 50ul #59472-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Periostin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	Periostin Antibody detects endogenous levels of total Periostin
Immunogen Description	A synthesized peptide derived from human Periostin
Other Names	Fasciclin-I like; OSF2; Periostin; Periostin osteoblast specific factor; PN; POSTN; RP11 412K4.1;
Accession No.	Uniprot:Q15063
Uniprot	Q15063
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 FC 1:50

## Images



Western blot analysis of Periostin expression in 293T cell lysate.

## Product Description

Binds to heparin. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization.

## Background

Binds to heparin. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization.

---

Note: This product is for in vitro research use only