

IP3 Receptor Rabbit mAb

Catalog No: #59496

Package Size: #59496-1 50ul #59496-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IP3 Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	IP3 Receptor Antibody detects endogenous levels of total IP3 Receptor
Immunogen Description	A synthesized peptide derived from human IP3 Receptor
Other Names	5-trisphosphate receptor; Inositol 1; InsP3R1; IP3; IP3 receptor; IP3R 1; IP3R; IP3R1; Itpr1; SCA15; SCA16; SCA29;
Accession No.	Uniprot:Q14643
Uniprot	Q14643
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of IP3 Receptor expression in HeLa cell lysate.

Product Description

Intracellular channel that mediates calcium release from the endoplasmic reticulum following stimulation by inositol 1,4,5-trisphosphate.

Background

Intracellular channel that mediates calcium release from the endoplasmic reticulum following stimulation by inositol 1,4,5-trisphosphate.

Note: This product is for in vitro research use only