

## RBP4 Rabbit mAb

Catalog No: #59508

Package Size: #59508-1 50ul #59508-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	RBP4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	RBP4 Antibody detects endogenous levels of total RBP4
Immunogen Description	A synthesized peptide derived from human RBP4
Other Names	PRBP; PRO2222; RBP; RBP4; RDCCAS;
Accession No.	Uniprot:P02753
Uniprot	P02753
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

## Images



Western blot analysis of RBP4 expression in HepG2 cell lysate.

## Product Description

Delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin, this prevents its loss by filtration through the kidney glomeruli.

## Background

---

Delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin, this prevents its loss by filtration through the kidney glomeruli.

---

Note: This product is for in vitro research use only