

TCTP Rabbit mAb

Catalog No: #59510

Package Size: #59510-1 50ul #59510-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCTP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	TCTP Antibody detects endogenous levels of total TCTP
Immunogen Description	A synthesized peptide derived from human TCTP
Other Names	Fortilin; Histamine releasing factor; HRF; Immunoglobulin E-dependent; p02; p23; TCTP; tpt1;
Accession No.	Uniprot:P13693
Uniprot	P13693
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:100

Images



Western blot analysis of TCTP Receptor expression in HeLa cell lysate.

Product Description

Involved in calcium binding and microtubule stabilization.

Background

Involved in calcium binding and microtubule stabilization.

Note: This product is for in vitro research use only