

## FOLR1 Rabbit mAb

Catalog No: #59515

Package Size: #59515-1 50ul #59515-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FOLR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	FOLR1 Antibody detects endogenous levels of total FOLR1
Immunogen Description	A synthesized peptide derived from human FOLR1
Other Names	Adult folate binding protein; FBP; Folate receptor 1; Folate Receptor 1 Precursor; Folate receptor alpha; FOLR1; FR alpha; KB cells FBP;
Accession No.	Uniprot:P15328
Uniprot	P15328
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of FOLR1 expression in HEK293 cell lysate.

## Product Description

Binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate to the interior of cells.

## Background

---

Binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate to the interior of cells.

---

Note: This product is for in vitro research use only