

GLUR Rabbit mAb

Catalog No: #59516

Package Size: #59516-1 50ul #59516-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GLUR Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	GLUR Antibody detects endogenous levels of total GLUR
Immunogen Description	A synthesized peptide derived from human GLUR
Other Names	GLR1; GLUR; Glutathione reductase; Glutathione reductase mitochondrial; Glutathione reductase, mitochondrial; GR; Gr1; GRase; GRD1; GSR;
Accession No.	Uniprot:P00390
Uniprot	P00390
Calculated MW	56kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of GLUR expression in (1) HeLa cell lysate; (2) Mouse brain lysate; (3) Rat kidney lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Maintains high levels of reduced glutathione in the cytosol.

Background

Maintains high levels of reduced glutathione in the cytosol.

Note: This product is for in vitro research use only